



# Unloader One® X

Clinically proven to relieve OA knee pain



Knee OA without bracing  
(bone-on-bone contact)



The 3-Point Leverage System  
(Unloader One® X shown)



Knee OA with bracing  
(space created between bones)

UNLOADER ONE® X

© ÖSSUR, 02. 2024, CRE-16913

## Peter Dudderidge's Unloader One® Brace Story



"My name is Peter Dudderidge, and I was referred to a Melbourne based Orthotics clinic for an Unloader One X knee brace due to painful knee osteoarthritis. From the initial consultation, where they took the time to understand my knee and leg issues, to diagnosing potential solutions and fitting the brace, the experience start to finish as been exceptional.

After a few months of using the Unloader One X knee brace in conjunction with heel-raising plates in my shoes, it has made a huge difference! The brace has made it possible for me to live, play cricket & sport in general is now a lot easier and bearable.

I have gone from thinking I needed an immediate knee replacement to feeling like a replacement is now a long way off."

- Peter Dudderidge

*Veterans cricketer Australia over 60's cricket captain 2023*

[Read Unloader One Stories](#)

# Unloader One<sup>®</sup> Upgrades

The Unloader brace range has undergone multiple upgrades to streamline the product experience and provide a consistent identity across the OA portfolio.

## Reduction in Waste

The Unloader One brace range now comes solely in black, reducing waste by streamlining color options and optimising resource efficiency.

## Reduction in Hinge Components

The Unloader One X now features a carbon hinge with just 4 components, offering a lightweight and robust solution. This represents a significant 70% reduction from 14 components, ensuring simplicity without sacrificing strength or durability. Additionally, tool-free flexion and extension stops eliminate the need for extra tools.

## New Sensil Technology

Sensil is printed onto thermoformed padding, molding precisely to the structure's fit and 3D shape. Silicone is printed on skin-contact areas to minimise waste.



Unloader One



Unloader One X

## Customer Care

Contact the Customer Care team today!

Email

Meet the team



Össur ANZ

**Sydney:** 26 Ross St, North Parramatta, NSW 2151 | 1300 123 268

**Melbourne:** 9 Market Dr, Bayswater North, VIC 3153 | 1300 123 268

**Auckland:** 6/80 Westpoint Dr, Hobsonville, Auckland, 0618 | 0800 369 524

[www.ossur.com.au](http://www.ossur.com.au) | [www.ossur.co.nz](http://www.ossur.co.nz)

© 2025 Össur

[Unsubscribe](#)

**Want to change the types of emails you receive from Össur?**

Selecting your preferences helps us tailor content to your interests, ensuring emails are relevant and engaging to you.

[Update your preferences](#)