



Hello,


The Unloader brace range has undergone multiple upgrades to streamline the product experience and provide a consistent identity across the OA portfolio. While all Unloader brace product codes and prices will remain the same, the product codes for accessories will progressively change as the range transitions to the new color. During this transition, you may also receive mixed inventory of old and new versions.

What are the upgrades?

1. Replace Sensil to Jersey Knit Thigh Liner
2. Removal Accessories & Polybag from product box
3. Replace Popliteal-, Hinge-, Ratchet Pads to Jersey Knit
4. Replace Current-strap material to UOX Strap Material
5. Tool-free Carbon Hinge & higher Quality Hinge Arms
6. Updated Sensil Calf Liner
7. Change from Grey to Black with Blue Details
8. Current Box is replaced with the Bracing Simplified Product Box

The decision to unify our product portfolio is based on our belief that consistency and uniformity provide our customers with an enhanced experience. By harmonising the colour scheme, we:

- Create an aesthetically appealing and cohesive appearance across our products.
- Minimise the need for various colour pigments and variants.
- Reduce waste and conserve resources



[Contact Us About Unloader](#)

Customer Care

Contact the Customer Care team today!

[Email](#)

[Meet the team](#)



Össur ANZ

Sydney: 26 Ross St, North Parramatta, NSW 2151 | 1300 123 268

Melbourne: 9 Market Dr, Bayswater North, VIC 3153 | 1300 123 268

Auckland: 6/80 Westpoint Dr, Hobsonville, Auckland, 0618 | 0800 369 524

www.ossur.com.au | www.ossur.co.nz

© 2025 Össur

[Unsubscribe](#)

Want to change the types of emails you receive from Össur?

Selecting your preferences helps us tailor content to your interests, ensuring emails are relevant and engaging to you.

[Update your preferences](#)