

Össur Leg

K-3 AK Order Form

Step 1 General Information

Prosthetist Name: _____ Acc Nr: _____
 Patient Identifier: _____ Weight: _____ kg Date of Birth: ____ / ____ / ____

Step 2 Select Suspension System



Seal-In®	Seal-In® Locking TF	Suspension Components	
ICEROSS SEAL-IN® X TF	ICEROSS SEAL-IN® X LOCKING TF	ICELOCK® CLUTCH 211	L-211000
Standard Conical	Standard Conical	OR ICELOCK® CLUTCH 4H 214	L-214000
Size: _____	Size: _____	OR *ICELOCK® 621 RATCHET	L-621000
Ring size: _____	Ring size: _____	OR ICELOCK® LANYARD 631	L-631000
Classic Volume Grip	Classic Volume Grip	OR ICELOCK® SMOOTH 651	L-651000
ICEROSS SEAL-IN® X5 TF	*ICEROSS® TF LOCKING LINER	OR *ICELOCK® EXPULSION VALVE 552	L-552000
Standard Conical	Standard Conical	OR UNITY TF VALVE KIT	UTF0001
Size: _____	Size: _____	OR ICELOCK 544 DIRECT SOCKET ADAPTER PLATE**	L-544050
*ICEROSS® TF SEAL-IN®		OR ICELOCK 544 PLATE EXPULSION**	L-544521
Standard Conical		OR ICELOCK 544 PLATE UNITY**	L-544511
Size: _____		OR ICELOCK 850 HYBRID FOR DIRECT SOCKET	L-850001
4 SEAL CLASSIC™		OR ICELOCK 850 HYBRID FOR 4-HOLE LAMINATION	L-850002
Standard Conical		OR ICELOCK 850 HYBRID FOR LAMINATED SOCKETS	L-850003
Size: _____			

Step 3 Select Socket - Direct Socket TF



Braid Selection

Liner Size	Braid Size
23.5 - 28	5"
30 - 32	5" OR 7" with Standard Distal Connector
34 - 40	7"*
45 - 55	9"

* Consider 7 inch braids with large distal connector during assessment of residual limb size, length and amount of soft tissue.

ICEROSS® CASTING LINER TF

Standard Conical

Size: _____

Direct Socket Material Kit Selection


5"	Glass Fibers with medium distal connector, Neutral Resin	M-501105
7"	Glass Fibers with medium distal connector, Neutral Resin	M-5001 07
7"	Glass Fibers with large distal connector, Neutral Resin	M-501107
9"	Glass Fibers with large distal connector, Neutral Resin	M-501109
5"	Glass Fibers with medium distal connector, Black Resin	M-501115
7"	Glass Fibers with medium distal connector, Black Resin	M-500117
7"	Glass Fibers with large distal connector, Black Resin	M-501117
9"	Glass Fibers with large distal connector, Black Resin	M-501119
5"	Basalt Fibers with medium distal connector, Black Resin	M-501315
7"	Basalt Fibers with medium distal connector, Black Resin	M-500317
7"	Basalt Fibers with large distal connector, Black Resin	M-501317
9"	Basalt Fibers with large distal connector, Black Resin	M-501319

Össur Brim TF Size Selection

Size	Circumference	
36	35 - 37 cm	M-FBS360
38	37 - 39 cm	M-FBS380
40	39 - 41 cm	M-FBS400
42	41 - 43 cm	M-FBS420
44	43 - 45 cm	M-FBS440
46	45 - 47.5 cm	M-FBS460
48.5	47.5 - 50 cm	M-FBS485
51	50 - 52.5 cm	M-FBS510
53.5	52.5 - 55 cm	M-FBS535
56	55 - 57.5 cm	M-FBS560
58.5	57.5 - 60 cm	M-FBS585
61	60 - 63 cm	M-FBS610
64	63 - 66 cm	M-FBS640
67	66 - 69 cm	M-FBS670
70	69 - 72 cm	M-FBS700



Fig. 1



BRIM SIZE SELECTION


1. Measure the circumference 1-2 cm below the level of ischial tuberosity (see fig.1)
2. Select the correct brim size from the table above.

Note. If the circumference is between two sizes, select the larger size. If the brim will be placed high on the lateral side, think of selecting one size larger.

Socket Components Selection

600 MODULAR SPACER FOR ICELOCK 600 SERIES	L-661000
BOTTOM SEAL FOR DIRECT SOCKET	M-400206

Step 4 Select Knee



*OHP3 KNEE™

*WITH MALE PYRAMID
WITH LOOP ADAPTER AND
FEMALE PYLON
WITH IKF ADAPTER

OP5 KNEE™

WITH MALE PYRAMID
WITH LOOP ADAPTER AND
FEMALE PYLON
WITH IKF ADAPTER

RHEO KNEE®

RHEO KNEE® XC
TOTAL KNEE® 2000
TOTAL KNEE® 2100
MAUCH® KNEE
CAPITAL KNEE
POWER KNEE™
NAVII® KNEE

ICON KNEE

PASO KNEE
WITH MALE PYRAMID
WITH LOOP ADAPTER AND
FEMALE PYLON
WITH IKF ADAPTER

Step 5 Select Foot + Cover



*VARI-FLEX®	▲PRO-FLEX® LP TORSION	▲PRO-FLEX® ST	*FOOT COVER	Beige	Brown
*▲VARI-FLEX® LP	▲PRO-FLEX® XC	▲PROPRIO FOOT®	Foot Size:		
▲PRO-FLEX® PIVOT	▲PRO-FLEX® XC TORSION	▲PRO-FLEX® TERRA	Left/Right:		
▲PRO-FLEX® LP	▲PRO-FLEX® LP ALIGN		Category:		

▲ Please note delivery period 14 working days for this product

Step 6 Select Adapter Combinations

SOCKET ADAPTERS

4-HOLE SOCKET ADAPTER AL 100kg	A-122100
4-HOLE SOCKET ADAPTER TI 166kg	A-125100
4-HOLE CONNECTION PLATE AL 136kg	1651014
4-PRONG SOCKET ADAPTER SS 166kg	A-114040
3-PRONG SOCKET ADAPTER SS 166kg	A-114030
MALE PYRAMID SOCKET ADAPTER TI 166kg	A-135100

PYRAMID ADAPTERS

PYRAMID INSERTS

MALE PYRAMID INSERT FOR PRONG TI 166kg	A-135300
FEMALE PYRAMID INSERT FOR PRONG TI 166kg	A-145300
FEMALE PYRAMID INSERT FOR PRONG Low Profile TI 166kg	A-145310
MALE PYRAMID INSERT FOR PRONG SS Insert 100kg	A-134300
FEMALE PYRAMID INSERT FOR PRONG SS 100kg	A-144300

4 HOLE ADAPTERS

4-HOLE MALE PYRAMID AL with SS Insert 100kg	A-233100
4-HOLE FEMALE PYRAMID AL 100kg	A-242100
4-HOLE TUBE CLAMP AL 100kg	A-322100
4-HOLE MALE PYRAMID WITH ROTATION SS 136kg	A-234320
4-HOLE FEMALE PYRAMID RECEIVER WITH ROTATION SS 136kg	A-244320
4-HOLE MALE PYRAMID TI 166kg	A-235100
4-HOLE MALE PYRAMID WITH ROTATION TI 166kg	A-235300
4-HOLE FEMALE PYRAMID TI 166kg	A-245100
4-HOLE FEMALE PYRAMID WITH ROTATION TI 166kg	A-245300
4-HOLE TUBE CLAMP TI 166kg	A-325100
EURO 4-HOLE ADAPTER AL 125kg	A-122300

SHIFT ADAPTERS

FEMALE PYRAMID TUBE CLAMP WITH 14MM OFFSET AL 100kg	A-342430
4-HOLE FEMALE PYRAMID WITH 20/20 OFFSET AL 125kg	A-242500
4-HOLE FEMALE PYRAMID 20/20 OFFSET 10° RIGHT AL 125kg	A-242510
4-HOLE FEMALE PYRAMID 20/20 OFFSET 10° LEFT AL 125kg	A-242520
MALE-FEMALE DOUBLE ADAPTER WITH 8MM OFFSET & 10° ANGLE SS 125kg	A-554700
4-HOLE MALE PYRAMID WITH 5/16" / 8MM OFFSET TI 166kg	A-235400
SHIFT ADAPTER WITH TITANIUM TI 166kg	A-645210
FLEXION PLATE 45MM 10°	A-282440
FLEXION PLATE 60MM 10°	A-282460
FLEXION PLATE 75MM 15°	A-282470

PYLONS AND TUBES

100kg

FEMALE PYLON KIT SHORT A-762111 AL 100kg	A-762111
FEMALE PYLON KIT LONG A-762121 AL 100kg	A-762121
HEIGHT ADJUSTABLE PYLON SHORT 77-95MM AL 100kg	A-542350
HEIGHT ADJUSTABLE PYLON LONG 94-128MM AL 100kg	A-542350
FEMALE PYLON SHORT AL 100kg	A-742010
FEMALE DOUBLE ADAPTER 32MM AL 100kg	A-542110
FEMALE DOUBLE ADAPTER 45MM AL 100kg	A-542120
FEMALE DOUBLE ADAPTER 60MM AL 100kg	A-542130
FEMALE DOUBLE ADAPTER 75 MM AL 100kg	A-542140
MALE DOUBLE ADAPTER AL with SS Pyramid Insert 100kg	A-533100
MALE PYRAMID TUBE CLAMP AL 100kg	A-333100
FEMALE PYRAMID TUBE CLAMP AL 100kg	A-342100
MALE SINGLE ADAPTER - SHORT AL 100kg	A-433110
MALE SINGLE ADAPTER - LONGAL 100kg	A-433120*
FEMALE SINGLE ADAPTER - SHORT AL 100kg	A-442100*
FEMALE SINGLE ADAPTER - LONG AL 100kg	A-442110*

166 KG

FEMALE PYLON SHORT TI 166kg	A-748030
FEMALE DOUBLE ADAPTER 32MM TI 166kg	A-545110
FEMALE DOUBLE ADAPTER 45MM TI 166kg	A-545120
FEMALE DOUBLE ADAPTER 60MM TI 166kg	A-545130
FEMALE DOUBLE ADAPTER 75 MM TI 166kg	A-545140
MALE DOUBLE ADAPTER TI 166kg	A-535100
MALE PYRAMID TUBE CLAMP TI 166kg	A-335100
FEMALE PYRAMID TUBE CLAMP TI 166kg	A-345100
MALE SINGLE ADAPTER - SHORT TI 166kg	A-435110
MALE SINGLE ADAPTER - LONG TI 166kg	A-435120
FEMALE SINGLE ADAPTER TI 166kg	A-445100

* Reimbursed for PMB