

TT Check Socket Order Form

Send all casts, sockets, files to: 9545 S John Young Parkway, Ste #100, Orlando, FL 32819
tel: (888)839-6213, fax: (800)788-9878 | cs@ossur.com | www.ossur.com

Billing & Shipping Information

Össur Account #:

Contact for ordering or delivery questions:

Contact:

Email:

Mobile Phone:

Ship To:

Company:

Address:

City/State/Zip:

Contact:

Phone:

Fax:

Email:

PO#:

Shipping Check Priority: Next Day Air **2nd Day**
Include return shipping label (\$50 charge inside 48 states)


Completion of this order form with the most accurate and up-to-date information, including all patient information, measurements, and construction information enables us to provide the highest quality product for you and your patient.

Required Patient Information

Last Name: Height: ft. in. Impact Level: High Med Low Amputation Level: TT Symes
First Name: Weight: lbs. K-Level: K1 K2 K3 K4 Left Right Bilat Order

Check Socket Measurements

Call to discuss before design

Measurements**			Shape	Socket Styles	Trim Lines
	Level	Measurement	Shape by cast Shape by app scan Shape by CAD file CAD file name:	Total Surface Bearing Patellar Tendon Bearing Modified PTB	Standard Supracondylar Supracondylar Suprapatellar
	10cm		Liner Used	*Angles Flexion: Adduction: *If a specific alignment isn't provided, ÖCS' standard alignment is 0° flexion and 0° adduction.	
	5cm				
	0cm		Total Reduction *3% global reduction (if using Iceross liners) OR 4% for any other liner unless requested differently. % *The percent reduction needed for correct socket fit is highly dependent on practitioner's measurement technique as well as patient's tissue density. It's not uncommon to require a different percent reduction. Ply		
	5cm				
	10cm				
	15cm				
	20cm				
	PML	cm			
	ML	cm			
AP	cm				
			Notes:		

**Length and circumferences are required. If not provided and socket fit requires remake, customer is responsible for cost of remake.

Cast sent to Össur Custom Solutions (optional) Check here to receive foam carving for this order.

TT Check Socket Order Form

Special Fabrications


Flexible Inner Socket Plastics	Fab Custom Inner From Sheet Material		Other
Flexible EVA with Silicone (Proflex with Silicone)	Lightweight Polyethylene – Firm 25 Shore A (PeLite)	Stovepipe	Lightweight Ethylene Vinyl Acetate (Keasy Cone) 6.5mm profile, 10mm distal
Flexible EVA without Silicone (Proflex without Silicone)	12mm (1/2") beige white	Slit in lateral side	
Flexible EVA – Black (Optek Flex Comfort - Black)	5mm (3/16") beige	Slit in medial side	
Flexible EVA – Optek Flex Comfort White + Black (BiLam)	3mm (1/8") beige		
Modified Polyethylene (MPE) EVA (ÖssurFlex)	Lightweight Polyethylene 25 Shore A (Bocklite) 6mm (1/4")		

Fab Custom Pad	RevoFit
<p>Lightweight Polyethylene – Medium Shore A (Plastizote)</p> <p>12mm (1/2") pink white</p> <p>6mm (1/4") pink</p> <p>5mm (3/16") pink</p> <p>Lightweight Polyethylene – Firm 25 Shore A (PeLite)</p> <p>12mm (1/2") beige white</p> <p>5mm (3/16") beige</p> <p>3mm (1/8") beige</p> <p>Location</p> <p>Distal end (1/2") Other (see notes below)</p>	<p>Dial Location</p> <p>Distal anterior medial Proximal lateral</p> <p>Other (see notes below)</p> <p>Panel Locations</p> <p>Pretibial & posterior (with pads)</p> <p>Posterior only (with pad)</p> <p>Pad on Panel</p> <p>Firm Ethylene Vinyl Acetate 55 Shore A (Puff)</p> <p>12mm (1/2") black</p> <p>10mm (3/8") black</p> <p>6mm (1/4") black</p>
Door	
Fab door with Velcro strap(s) (draw door on cast/check socket/use notes below)	

Notes

TT Check Socket Order Form

Suspension & Socket Adapter					
Locking Options					
Icelock 621 Ratchet & Extra Pin  	Icelock 651 Smooth & Extra Pin  	Icelock 562 Hybrid  	Icelock 600 XM 	Icelock 211 Clutch & Extra Pin Icelock 214 4-Hole Clutch & Extra Pin 	Icelock 234 Lanyard 
Works With					
Icelock 672 Male Pyramid Ti  Icelock 673 Male Pyramid Al 	Icelock 674 4-Hole Adapter Al Icelock 675 4-Hole Adapter Ti  				
Expulsion Options			Unity Elevated Vacuum Options		
Icelock 544 Expulsion Kit 	Icelock 551 Expulsion Valve  	Unity Valve TT  		Icelock 544 Unity Plate 	
Icelock 544 Expulsion Valve Kit for sleeves 	Works With				
	Icelock 544 Socket Adapter 	3-Prong Socket Adapter 4-Prong Socket Adapter			
Other (write in):					

 Waterproof rated Össur products provide permanent protection against harmful ingress of fresh, salt or chlorinated water, even after submersion.

TT Definitive Socket Order Form

Send all casts, sockets, files to: 9545 S John Young Parkway, Ste #100, Orlando, FL 32819
tel: (888)839-6213, fax: (800)788-9878 | cs@ossur.com | www.ossur.com

Billing & Shipping Information

Ossur Account #:

CONTACT FOR ORDERING OR DELIVERY QUESTIONS:

Contact:

Ship To:

Company:

Address:

City/State/Zip:

Contact:

Phone:

Email:

PO#:

Shipping Check Priority:

Next Day Air

2nd Day

Completion of this order form with the most accurate and up-to-date information, including all patient information, measurements, and construction information enables us to provide the highest quality product for you and your patient.

Patient Information

Patient Name/Reference:

Amputation Level: BK SYMES

Activity Level: 1 2 3 4


Weight:

Left: Right: Bilat Order:

Definitive Socket

Call to discuss before fabrication

TT Measurement Chart



From Socket

From Cast

From CAD/AOP file dated:

Lamination Finish

Carbon

Basalt

Color (PRS Pigment 1-16)

Decorative Lamination

Alignment

Transfer

Transfer and zero out

Use Lines

Other

Use Angles

Flexion

Adduction

Modifications

Socket Sizing

Prosthetic Height

Increase
Decrease

%:

Ply:

(optional)

Transfer

Extend:

Shorten:

TT Definitive Socket Order Form

Special Fabrications

Flexible Inner Socket Plastics	Fab Custom Inner From Sheet Material		Other
Flexible EVA with Silicone (Proflex with Silicone)	Lightweight Polyethylene – Firm 25 Shore A (PeLite)	Stovepipe	Lightweight Ethylene Vinyl Acetate (Keasy Cone) 6.5mm profile, 10mm distal
Flexible EVA without Silicone (Proflex without Silicone)	12mm (1/2") beige white	Slit in lateral side	
Flexible EVA – Black (Optek Flex Comfort - Black)	5mm (3/16") beige	Slit in medial side	
Flexible EVA – Optek Flex Comfort White + Black (BiLam)	3mm (1/8") beige		
Modified Polyethylene (MPE) EVA (ÖssurFlex)	Lightweight Polyethylene 25 Shore A (Bocklite) 6mm (1/4")		

Fab Custom Pad	RevoFit
<p>Lightweight Polyethylene – Medium Shore A (Plastizote)</p> <p>12mm (1/2") pink white</p> <p>6mm (1/4") pink</p> <p>5mm (3/16") pink</p> <p>Lightweight Polyethylene – Firm 25 Shore A (PeLite)</p> <p>12mm (1/2") beige white</p> <p>5mm (3/16") beige</p> <p>3mm (1/8") beige</p> <p>Location</p> <p>Distal end (1/2") Other (see notes below)</p>	<p>Dial Location</p> <p>Distal anterior medial Proximal lateral</p> <p>Other (see notes below)</p> <p>Panel Locations</p> <p>Pretibial & posterior (with pads)</p> <p>Posterior only (with pad)</p> <p>Pad on Panel</p> <p>Firm Ethylene Vinyl Acetate 55 Shore A (Puff)</p> <p>12mm (1/2") black</p> <p>10mm (3/8") black</p> <p>6mm (1/4") black</p> <p>See markings on cast/socket</p> <p>Other (see notes below)</p>
<p>Custom shaped foam cover (provide meas. in diagram on previous page)</p> <p>Fab cosmetic distal end</p> <p>Fab windows, cutouts or Dacron strap (draw windows or cutouts on check socket/use notes below)</p> <p>Fab door with Velcro strap(s) (draw door on check socket/use notes below)</p>	

Notes