



Mechanical Upper Limb Scanning Now Available!



You Shared, We Listened

Mechanical Upper Limb scanning is here!

Following the recent brand name update from Naked Prosthetics to Össur Upper Limb, scanning is now available for Mechanical Upper Limb devices.

We have listened to your feedback and made improvements to streamline intake, reduce delays, and help lower the risk of impression rejection.

Scanning Guide

What has changed

- Digital hand scans can now be submitted for Mechanical Upper Limb devices
 - Accepted file formats: PLY, STL, OBJ, MF
- Create scans using conventional scanners already in your clinic
 - Plaster molds remain available if preferred

What this means in practice

- Less time spent preparing and shipping molds
- Reduced transit delays, especially for international shipments
 - Lower risk of impression rejection
- A more efficient intake process from start to finish

Scanning Intake Made Simple

For practical information on scanning for MCPDriver and ThumbDriver, view this short course to help you get started with confidence.

MCPDriver & ThumbDriver Scanning Intake Course

New Photo Intake Mat

A new Photo Intake Mat is now available to support clearer, more consistent clinical images when ordering Mechanical Upper Limb devices.

The mat helps achieve correct scaling and image clarity, reducing retakes and supporting a smoother ordering process.


- Item number NPD1 0144
- Available at no cost

**Request yours through your
Clinical Specialist.**

[Photo Intake Mat User Guide](#)

Our clinical team will provide in field guidance and support to help you adopt digital scanning with confidence.

To arrange a setup session or ask a question, please contact your local Össur clinical specialist.



Clinical Team

Contact the Clinical Team today!

Email

Meet the team



Össur ANZ

Sydney: 26 Ross St, North Parramatta, NSW 2151 | 1300 123 268

Melbourne: 9 Market Dr, Bayswater North, VIC 3153 | 1300 123 268

Auckland: 6/80 Westpoint Dr, Hobsonville, Auckland, 0618 | 0800 369 524

www.ossur.com.au | www.ossur.co.nz

© 2026 Össur

[Unsubscribe](#)

Want to change the types of emails you receive from Össur?

Selecting your preferences helps us tailor content to your interests, ensuring emails are relevant and engaging to you.

[Update your preferences](#)