

Össur Leg

K-1 BK Order Form

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Step 1 General Information

Prosthetist Name: _____ Acc Nr: _____
 Patient Identifier: _____ Weight: _____ kg Date of Birth: ____ / ____ / _____

Step 2 Select Suspension System



| Seal-In® | Cushion | Locking | PLUS Suspension System |
|--|--|---|--|
| ICEROSS SEAL-IN® X 3mm 6mm Size: _____ Ring size: _____ Classic Volume Grip | ICEROSS® SYNERGY™ CUSHION Wave Size: _____ ICEROSS DERMO® CUSHION 3mm 6mm 3mm Wave Size: _____ | ICEROSS® SYNERGY™ LOCKING Wave Size: _____ ICEROSS DERMO® LOCKING 3mm 6mm 3mm Wave 3mm Conical Wave Size: _____ | ICELOCK ICELOCK® CLUTCH 211 L-211000 OR ICELOCK® CLUTCH 4H 214 L-214000 OR *ICELOCK® 621 RATCHET L-621000 OR ICELOCK® LANYARD 631 L-631000 OR ICELOCK® SMOOTH 651 L-651000 OR ICELOCK® 562 HYBRID L-562000 OR ICELOCK® 562 HYBRID FOR PASSIVE VACUUM L-562500 OR ICELOCK® 562 HYBRID DS L-562400 OR ICELOCK® 562 HYBRID FOR PASSIVE VACUUM DS L-562410 OR ICELOCK® 850 HYBRID DS L-850001 OR ICELOCK® 850 HYBRID 4H L-850002 OR ICELOCK® 850 HYBRID LS L-850003 |
| ICEROSS SEAL-IN® V Standard Profile High Profile Size: _____ | *ICEROSS COMFORT® CUSHION 3mm 6mm Size: _____ | *ICEROSS COMFORT® LOCKING 3mm 6mm 3mm Wave Size: _____ | PLUS Expulsion VALVE SYSTEM *ICELOCK® EXPULSION VALVE 551 L-551002 |
| ICEROSS SEAL-IN® X5 Wave Size: _____ *ICEROSS SEAL-IN DERMO® Size: _____ | PLUS SLEEVE ICEROSS® SLEEVE OR ICEFLEX® BALANCE OR ICEFORM® SLEEVE OR *GENU SLEEVE Sleeve size: _____ | | OR PLUS UNITY SYSTEM UNITY TT VALVE KIT UTT0001 UNITY PUMP |
| ICEROSS SEAL-IN X LOCKING 3mm 6mm Classic Volume Grip Ring size: _____ | | | |

Step 3 Select Socket - Direct Socket TT



Braid Selection

| Liner Size | Braid Size |
|------------|------------|
| 16 - 22 | 4" |
| 23.5 - 30 | 5" |
| 32 - 34 | 7" |

ICEROSS® CASTING LINER TT

Size: _____

Socket Components Selection

| | |
|---|----------|
| 600 MODULAR SPACER FOR ICELOCK 600 SERIES | L-661000 |
| BOTTOM SEAL FOR DIRECT SOCKET | M-400206 |

Direct Socket Material Kit Selection

| | |
|---|----------|
| 4" Glass Fibers with small distal connector, Neutral Resin | M-500104 |
| 5" Glass Fibers with medium distal connector, Neutral Resin | M-500105 |
| 7" Glass Fibers with medium distal connector, Neutral Resin | M-500107 |
| 7" Glass Fibers with large distal connector, Neutral Resin | M-501107 |
| 4" Glass Fibers with small distal connector, Black Resin | M-500114 |
| 5" Glass Fibers with medium distal connector, Black Resin | M-500115 |
| 7" Glass Fibers with medium distal connector, Black Resin | M-500117 |
| 7" Glass Fibers with large distal connector, Black Resin | M-501117 |
| 5" Carbon Fiber with medium distal connector, Black Resin | M-500205 |
| 7" Carbon Fiber with medium distal connector, Black Resin | M-500207 |
| 4" Basalt Fibers with small distal connector, Black Resin | M-500314 |
| 5" Basalt Fibers with medium distal connector, Black Resin | M-500315 |
| 7" Basalt Fibers with medium distal connector, Black Resin | M-500317 |
| 7" Basalt Fibers with large distal connector, Black Resin | M-501317 |

Step 4 Select Foot + Cover



*FLEX-FOOT® BALANCE



BALANCE™ FOOT S



BALANCE™ FOOT J



BREEZE™ FOOT



*FOOT COVER

Beige

Brown

Foot Size: _____

Left/Right: _____

Category: _____

Step 5 Select Adapter Combinations

SOCKET ADAPTERS

| | |
|--------------------------------------|----------|
| 4-HOLE SOCKET ADAPTER AL 100kg | A-122100 |
| 4-HOLE SOCKET ADAPTER TI 166kg | A-125100 |
| 4-HOLE CONNECTION PLATE AL 136kg | 1651014 |
| 4-PRONG SOCKET ADAPTER SS 166kg | A-114040 |
| 3-PRONG SOCKET ADAPTER SS 166kg | A-114030 |
| MALE PYRAMID SOCKET ADAPTER TI 166kg | A-135100 |

PYRAMID ADAPTERS

PYRAMID INSERTS

| | |
|--|----------|
| MALE PYRAMID INSERT FOR PRONG TI 166kg | A-135300 |
| FEMALE PYRAMID INSERT FOR PRONG TI 166kg | A-145300 |
| FEMALE PYRAMID INSERT FOR PRONG Low Profile TI 166kg | A-145310 |
| MALE PYRAMID INSERT FOR PRONG SS Insert 100kg | A-134300 |
| FEMALE PYRAMID INSERT FOR PRONG SS 100kg | A-144300 |

4 HOLE ADAPTERS

| | |
|---|----------|
| 4-HOLE MALE PYRAMID AL with SS Insert 100kg | A-233100 |
| 4-HOLE FEMALE PYRAMID AL 100kg | A-242100 |
| 4-HOLE TUBE CLAMP AL 100kg | A-322100 |
| 4-HOLE MALE PYRAMID WITH ROTATION SS 136kg | A-234320 |
| 4-HOLE FEMALE PYRAMID RECEIVER WITH ROTATION SS 136kg | A-244320 |
| 4-HOLE MALE PYRAMID TI 166kg | A-235100 |
| 4-HOLE MALE PYRAMID WITH ROTATION TI 166kg | A-235300 |
| 4-HOLE FEMALE PYRAMID TI 166kg | A-245100 |
| 4-HOLE FEMALE PYRAMID WITH ROTATION TI 166kg | A-245300 |
| 4-HOLE TUBE CLAMP TI 166kg | A-325100 |
| EURO 4-HOLE ADAPTER AL 125kg | A-122300 |

SHIFT ADAPTERS

| | |
|---|----------|
| FEMALE PYRAMID TUBE CLAMP WITH 14MM OFFSET AL 100kg | A-342430 |
| 4-HOLE FEMALE PYRAMID WITH 20/20 OFFSET AL 125kg | A-242500 |
| 4-HOLE FEMALE PYRAMID 20/20 OFFSET 10° RIGHT AL 125kg | A-242510 |
| 4-HOLE FEMALE PYRAMID 20/20 OFFSET 10° LEFT AL 125kg | A-242520 |
| MALE-FEMALE DOUBLE ADAPTER WITH 8MM OFFSET & 10° ANGLE SS 125kg | A-554700 |
| 4-HOLE MALE PYRAMID WITH 5/16" / 8MM OFFSET TI 166kg | A-235400 |
| SHIFT ADAPTER WITH TITANIUM TI 166kg | A-645210 |

PYLONS AND TUBES

100kg

| | |
|--|----------|
| FEMALE PYLON KIT SHORT A-762111 AL 100kg | A-762111 |
| FEMALE PYLON KIT LONG A-762121 AL 100kg | A-762121 |

| | |
|--|----------|
| HEIGHT ADJUSTABLE PYLON SHORT 77-95MM AL 100kg | A-542350 |
| HEIGHT ADJUSTABLE PYLON LONG 94-128MM AL 100kg | A-542350 |

| | |
|--------------------------------------|----------|
| FEMALE PYLON SHORT AL 100kg | A-742101 |
| FEMALE DOUBLE ADAPTER 32MM AL 100kg | A-542110 |
| FEMALE DOUBLE ADAPTER 45MM AL 100kg | A-542120 |
| FEMALE DOUBLE ADAPTER 60MM AL 100kg | A-542130 |
| FEMALE DOUBLE ADAPTER 75 MM AL 100kg | A-542140 |

| | |
|---|----------|
| MALE DOUBLE ADAPTER AL with SS Pyramid Insert 100kg | A-533100 |
|---|----------|

| | |
|------------------------------------|----------|
| MALE PYRAMID TUBE CLAMP AL 100kg | A-333100 |
| FEMALE PYRAMID TUBE CLAMP AL 100kg | A-342100 |

| | |
|--|-----------|
| MALE SINGLE ADAPTER - SHORT AL 100kg | A-433110 |
| MALE SINGLE ADAPTER - LONGAL 100kg | A-433120* |
| FEMALE SINGLE ADAPTER - SHORT AL 100kg | A-442100* |
| FEMALE SINGLE ADAPTER - LONG AL 100kg | A-442110* |

166 KG

| | |
|--------------------------------------|----------|
| FEMALE PYLON SHORT TI 166kg | A-748030 |
| FEMALE DOUBLE ADAPTER 32MM TI 166kg | A-545110 |
| FEMALE DOUBLE ADAPTER 45MM TI 166kg | A-545120 |
| FEMALE DOUBLE ADAPTER 60MM TI 166kg | A-545130 |
| FEMALE DOUBLE ADAPTER 75 MM TI 166kg | A-545140 |

| | |
|------------------------------|----------|
| MALE DOUBLE ADAPTER TI 166kg | A-535100 |
|------------------------------|----------|

| | |
|------------------------------------|----------|
| MALE PYRAMID TUBE CLAMP TI 166kg | A-335100 |
| FEMALE PYRAMID TUBE CLAMP TI 166kg | A-345100 |

| | |
|--------------------------------------|----------|
| MALE SINGLE ADAPTER - SHORT TI 166kg | A-435110 |
| MALE SINGLE ADAPTER - LONG TI 166kg | A-435120 |
| FEMALE SINGLE ADAPTER TI 166kg | A-445100 |

* Reimbursed for PMB



WWW.OSSUR.COM/EN-ZA

Össur South Africa
Unit 1, Stonewood
Business Park
9 Tanzanite Crescent
Kraaifontein
Cape Town
7550

T 0800 996 220
F 0800 981 744
E ossursa@ossur.com

