











ÖSSUR LEGS TRANSFEMORAL PROSTHESIS ORDER FORM


FACILITY INFORMATION	
Company Name: <input type="text"/>	Ship to Address: <input type="text"/>
Clinician Name: <input type="text"/>	
Contact Number: <input type="text"/>	
Email: <input type="text"/>	
Purchase Order# <input type="text"/>	

REQUIRED PATIENT INFORMATION	
Patient Name: <input type="text"/>	Impact Level: High Med Low
	Weight: kg. K-Level: K1 K2 K3 K4 Left Right

1. TF LINERS (check one then fill in desired sizes and options)
Limit 2 liners of same type, 2 Seal-Rings, 2 Distal Cups, 1 Sleeve - per order.

LINER OPTIONS									
TF Liners									
Iceross Casting Liner*	Iceross TF Liner	AKOS	Relax TF	Relax 3C TT Cushion	4Seal Classic	Seal-In X5	Seal-In HSM	Seal-In X	Seal-In X TF Locking
									
* casting liner is required for direct socket casting									
















LINER SPECIFICATIONS			
Liner Size		Optional: Add additional liner of same type	
<input type="text" value="LINER SIZE"/>	Standard Conical	<input type="text" value="LINER SIZE"/>	Standard Conical

OPTIONAL ADD-ONS			
Optional: Add Distal Cups (Limit 2)		Optional: Add Seal-Rings (Limit 2)	
<input type="text" value="DISTAL CUP SIZE"/>	Distal Cup 1	Seal-Ring 1	Seal-Ring 2
<input type="text" value="DISTAL CUP SIZE"/>	Distal Cup 2	 <input type="text" value="RING SIZE"/>	 <input type="text" value="RING SIZE"/>

NOTES:














2. SUSPENSION OPTIONS (Choose for conventional socket fabrication) Refer to section 3 for Direct Socket Suspension Options
Limit 1 lock/valve kit per order.

LOCKING OPTIONS

680 Fabrication Adaption Kit 	621 Ratchet Pin  	651 Smooth Pin  	214 Clutch  	211 Clutch  	200 Series Lanyard  
Icelock 850 Hybrid for 4-Hole Socket Add Optional Value: 551 Valve 552 Valve Unity TF Valve Kit  			Icelock 850 Hybrid for Laminated Socket Add Optional Value: 551 Valve 552 Valve Unity TF Valve Kit  		











Add Extra Pin: Long Standard Short 850 Hybrid Lock Pin

EXPULSION OPTIONS

544 Expulsion Plate Kit for Seal-In 	544 Expulsion Valve Kit for Sleeves 	Expulsion Valve 551 Expulsion Valve 552  	UNITY Unity TT Valve Kit Unity TF Valve Kit  	544 Expulsion Plate Kit with Unity   Unity TF Tube Kit  
		Works With		
		544 Socket Adaptor 	3-Prong 4-Prong  	

3. DIRECT SOCKET SUSPENSION OPTIONS (Choose when fabricating your socket using a Direct Socket Kit)
Limit 1 lock/valve kit per order








LOCKING

621 Ratchet Pin   Requires: Modular Spacer Plate 	651 Smooth Pin   Requires: Modular Spacer Plate 	214 Clutch Lock  	850 Hybrid Lock for Direct Socket   Add Optional Value: 551 Valve 552 Valve Unity TF Valve Kit
---	--	--	---

Add Extra Pin: Long Standard Short 850 Hybrid Lock Pin

3. DIRECT SOCKET SUSPENSION OPTIONS (Choose when fabricating your socket using a Direct Socket Kit)

Limit 1 lock/valve kit per order

EXPULSION		UNITY	
Expulsion Valve 551 Expulsion Valve 552  	Expulsion Plate 	Unity TT Valve Unity TF Valve  	544 Unity Plate  Unity TF Tube Kit 
Requires: Bottom Seal 	Requires: DS Adapter Plate 	Requires: Bottom Seal 	Requires: DS Adapter Plate 


4. OPTIONAL: DIRECT SOCKET TF KITS

Limit 1 Direct Socket Kit per order

GLASS				BASALT			
 		DISTAL CONNECTOR  Standard  Large		 		DISTAL CONNECTOR  Standard  Large	
Resin:	Black	Neutral		Black Resin:	5 inch	7 inch	9 inch
Braid:	5 inch	7 inch	9 inch				

BRIM (Select one when using a Direct Socket TF Kit)

Limit 1 with purchase of Direct Socket TF Kit

	Size:	36	38	45	40	42	44	46	48.5
		51	53.5	56	58.5	61	64	67	70

5. OPTIONAL: CONNECT TF PACKAGE KIT

Limit 1 Connect TF Socket Kit per order

Connect TF Kit

Select your Size:

- Medium Short
- Medium Standard
- Large Short
- Large Standard

Left Right





















































































Note: Connect TF is only available as an option when purchased with Icross TF Locking liner(s), AKOS, and Relax. One Connect TF socket kit per order. Connect TF includes integrated 200 series clutch lock.

6. ENDOSKELETAL COMPONENTS (check all that apply)

Limits apply on some components. Require a minimum of 2 units; Max 5 total components. Refer to below items;

Socket components: Limit 2; Compatible types
Pylons: Limit 2




























































Tube Clamps: Quantity limitations vary depending on choice of pylons or pylon kits
Torsion Shock Adapter: Limit 1 (when foot does not include torsion/shock adapter)

4-HOLE ADAPTERS		DOUBLE ADAPTERS		SINGLE ADAPTERS	
	4-Hole Female  		Female Double Adapter 32mm  		Male Single Adapter Short  
	4-Hole Flexion Plate 45mm offset 10° Flexion 		Female Double Adapter 45mm  		Male Single Adapter Long  
	4-Hole Flexion Plate 60mm offset 10° Flexion 		Female Double Adapter 60mm  		Female Single Adapter Short  
	4-Hole Flexion Plate 75mm offset 15° Flexion 		Female Double Adapter 75mm  		Female Single Adapter Long 
	4-Hole Female w/ Rotation  	TORSION SHOCK ADAPTER		TUBE CLAMPS	
	4-Hole Female Pyramid 20/20 Offset 10° Right 		Cat 1 		Male Tube Clamp  
	4-Hole Female Pyramid 20/20 Offset 10° Left 		Cat 2 		Female Tube Clamp  
			Cat 3 		4-Hole Tube Clamp  
PRONG ADAPTERS		PYLONS		HEIGHT ADJUSTABLE PYLONS	
	4-Prong Socket Adapter 		Female Short  		Standard 
	3-Prong Socket Adapter 		Female Kit Short  		Long 
	Male Pyramid Insert for Prong  		Female Long  	<p>Please tick this box if you would prefer Titanium adaptors.</p> <p>* Patients under 100kg will be supplied with Aluminium adaptors unless specified.</p> <p>* Patients between 100-166kg will be provided with Titanium adaptors.</p> <p>Please tick this box if you would prefer waterproof version of components, availability indicated by icons</p>	
	Female Pyramid Insert for Prong  		Female Kit Long  		
	Female Pyramid Insert for Prong LP 		Waterproof HD Titanium version of the above		
SPECIFY OTHER COMPONENTS AND/OR QUANTITIES DESIRED					

NOTES:

























































7. KNEE (check one and fill out the options below)

Limit 1 per order. Trial units are available for use prior to purchase. Bionic and Mechanical Knee units may not be returned or exchanged following purchase.

LOW ACTIVITY	MODERATE ACTIVITY	HIGH ACTIVITY
 <p>Locking Knee </p>	 <p>OP4 Knee™ </p>	 <p>Cheetah® Knee </p>
 <p>Balance™ Knee OFM1  Pyramid adapter  IKF adapter  Loop Adapter </p>	 <p>OP5 Knee™  Pyramid adapter  IKF adapter  Loop Adapter </p>	 <p>Total Knee® 2100  4-Prong Socket Adapter  Euro 4-Hole Adapter  Male PYR ADP (100kg)  Male PYR ADP (166kg)  Female PYD ADP </p>
 <p>Balance™ Knee OFM2 </p>	 <p>OHP3 Knee™  Pyramid adapter  IKF adapter  Loop Adapter </p>	
 <p>Total Knee® 1900  4-Prong Socket Adapter  Euro 4-Hole Adapter  Male PYR ADP (100kg)  Male PYR ADP (166kg)  Female PYR ADP </p>	 <p>Mauch® Knee (136kg) </p>	
	 <p>Total Knee® 2000  4-Prong Socket Adapter  Euro 4-Hole Adapter  Male PYR ADP (100kg)  Male PYR ADP (166kg)  Female PYR ADP </p>	<p style="text-align: center;">BIONIC SOLUTIONS</p> <p style="text-align: center;">MODERATE ACTIVITY</p>  <p>RHEO KNEE® </p>
LOW TO MODERATE ACTIVITY	LOW TO HIGH ACTIVITY	HIGH ACTIVITY
 <p>Power Knee™ </p>	 <p>Paso Knee™  Pyramid adapter  IKF adapter  Loop Adapter </p>	 <p>RHEO KNEE® XC </p>
NOTES:		
		 <p>NAVII® </p>

8. FOOT: WITH OR WITHOUT UNITY (check one and fill out the options below)

Limit 1 per order. Trial Proprio Foot units are available for use prior to purchase. Bionic units may not be returned or exchanged following purchase. Mechanical feet may be exchanged during the 60-day trial period.

FOOT OPTIONS:		Foot Size:	cm.	Category (1-9)	Left	Right	Foot Cover:	Beige	Brown
LOW ACTIVITY			MODERATE ACTIVITY			HIGH ACTIVITY			
	Flex-Foot® Balance with Unity®				LP Vari-Flex® with Unity®	 		Cheetah® Xplore Foot Cover Kit Lamination Connector Kit	
	Flex-Foot® Balance with D/P Flexion with Unity®				Vari-Flex®** with Unity®	 		Re-flex Shock™ with Unity® Standard Tall	 
	Balance™ Foot J with Unity®	 			Pro-Flex® Pivot with Unity®	 		Re-flex Rotate™ with Unity® Standard Male Tall Male Standard Female Tall Female	 
	Flex-Foot® Assure** with Unity®	 	*		Pro-Flex® LP with Unity®	 		Pro-Flex® XC Torsion with Unity®	 
	Balance™ Foot S with Unity®	 			Pro-Flex® LP Align			Pro-Flex® XC with Unity®	 
	Balance™ Foot S Torsion with Unity®	 		Other (Please Specify)				Pro-Flex® LP Torsion with Unity®	 
				<div style="border: 1px solid gray; height: 100px; width: 100%;"></div>				Pro-Flex® Terra with Unity®	 
BIONIC SOLUTION									
LOW TO MODERATE ACTIVITY			LOW TO HIGH ACTIVITY						
	PROPRIO FOOT® with Unity®	 			Pro-Flex® ST with Unity®	 		Pro-Flex® Modular	

* Waterproof only when used with tube kit

FOOT AND ADAPTER CHARACTERISTICS



WEATHERPROOF

A weatherproof device allows use in a wet and/or humid environment but does not allow submersion. Fresh water splashing against the enclosure from any direction shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.



WATERPROOF

Whether it's fresh water, salt water or chlorinated water, Waterproof Össur products will give you the maximum protection you need. The use of first-rate materials provides permanent protection against harmful ingress of water, even after submersion, whilst continuing to provide the functionality you expect.

NOTES: