



# i-Digits<sup>®</sup> Quantum



## The Future of Dexterity

**i-Digits Quantum is an electronic prosthesis that moves with you. It uses sensors on your forearm muscles to open and close the fingers as needed for daily tasks, whether you're using your phone, preparing a meal, or making a coffee.**

Designed for people missing one to five fingers, i-Digits is custom-fitted to your hand and works together with any remaining fingers.

### Features That Make Life Easier

- **Gesture Control** - A convenient way to switch grips. Simply move your i-Digits Quantum in one of four directions to quickly select the grip you need. This gives you fast access to grips especially useful when you're busy or on the go.
- **Conformable Grip** - The fingers move until they meet the object you're picking up, allowing them to naturally shape around it. This helps you grasp items of different sizes and forms, from a water bottle to a shopping bag.
- **Speed Boost** - The speed of the finger movement can be increased by up to 30%, allowing you to adjust it to best meet your needs.



Intelligent Motion  
Gesture Control



Triggers  
Muscle Control



Quick Grips<sup>®</sup>  
App Control



Grip Chips<sup>®</sup>  
Proximity Control