

# MCPDriver®

Restore Your Grip. Reclaim Your Routine.



**Designed for those missing a finger close to the base, the MCPDriver moves with your natural motion so you can get back to tasks that matter.**

### Your Movement is the Power

The MCPDriver is powered by you. By using the movement of your remaining finger, the device mimics the way an intact finger moves. When you move, it moves, giving you simple, intuitive control.

### A Fit as Unique as Your Hand

We don't believe in "one size fits all" - precision built to your exact measurements.

Built for Comfort - designed to be comfortable enough to wear from your morning coffee to your evening walk.

### Return to two handed activities

The MCPDriver makes it easier to rely on both hands again:

- **Unshakeable Strength:** Whether it's a power tool or a heavy suitcase, the MCPDriver provides the stability you need for demanding work.
- **Strong & Gentle Grip:** Get a firm, reliable hold on a steering wheel, a heavy tool, or a delicate glass.




- **Stay Connected:** An optional **touchscreen tip** can be added to your device, to allow you to use your smartphone or tablet.

### Ready for the Real World


Your prosthesis should be rugged enough for your toughest projects, not the other way around.

- **Waterproof & Durable:** Whether you're washing dishes, gardening in the rain, or building in a workshop, the MCPDriver is built to get wet and get dirty.
- **Personalised Style:** Choose from a wide range of colours and textures to make the device truly yours.

### Colours



### Textures



### Fasteners



Please note that due to the manufacturing process, variations in the colour, shade, and tone of nylon components and anodised fasteners may occur.