

PIPDriver®

Protect Your Finger. Regain Your Grip.



Small but Mighty. The PIPDriver is designed to protect your sensitive fingertip while giving you back the finesse needed to type, pick up a coin, or button a shirt.

Your Movement, Your Control

The **PIPDriver** works by following the bending and straightening of your remaining finger joint. When you move, it moves, giving you simple, intuitive control.

Made Just for You

Precision-built for your exact measurements, giving you a comfortable fit all day long.

Designed to Protect While Helping You Do More

Many people with partial finger loss experience sensitivity in the remaining part of the finger. The **PIPDriver** is designed to protect your finger while giving you back your grip. It gently shields the sensitive area from bumps, impact, and abrasion, helping you return to everyday activities.

Precision for the Smallest Details

The **PIPDriver** makes it easier to rely on both hands again, whether you're cooking, typing, opening packages, or carrying groceries.



It supports multiple grips and pinches, helping you handle everyday objects naturally, from a pen to a garden tool.

Easy to Live With


Your prosthesis should fit into your life, not the other way around.

- Easy to clean and care for
- Simple to put on and take off
- Strong, durable materials built with advanced 3D printing technology
- Designed to withstand fresh, salt, and chlorinated water - ready for real-life environments, not just indoor use
- Optional **touchscreen tip** lets you use your smartphone, tablet, or ATM

Colours



Textures



Fasteners



Please note that due to the manufacturing process, variations in the colour, shade, and tone of nylon components and anodised fasteners may occur.